

AccuPower[™] Ligation PreMix

AccuPower[™] Ligation PreMix is a new, powerful, ready-to-use ligation reagent optimized for more accurate ligation reactions.

- Convenient lyophilized ligation premix which consists of T4 DNA ligase, ATP, and reaction buffer
- Requires just 5 minutes for cohesive end DNA ligation
- Ligation activity and stability up to 4 months at room temperature

DESCRIPTION

The *AccuPower*^M Ligation PreMix is an all-in-one tube containing a lyophilized premix of all the necessary components for a ligation reaction. The mix includes T4 DNA ligase, ATP, reaction buffer and stabilizer. The chemical stabilizer (patent pending) maintains activity for 4 months even when stored at room temperature (25°C) and 2 years at -20°C. The *AccuPower*^M Ligation PreMix is designed for cohesive end ligation in 5 minutes and blunt end ligation in 10 minutes at room temperature. Rapid ligation is based on the combination of T4 DNA Ligase with an optimized premix buffer.

FEATURES AND BENEFITS

High Speed	The ligation can be performed in 5 minutes for cohesive ends and in 10 minutes for blunt end DNA at room temperature.	
Simplicity	The <i>AccuPower™</i> Ligation PreMix is ready to use, thus offering virtually no reaction set up.	
Stability	The chemical stabilizer maintains enzyme activity for 4 months at room temperature, and 2 years at -20°C.	
Flexibility	Suitable for all common ligation reactions.	
Reproducibility	Batch manufacturing under strict ISO 9001 quality control system.	

AccuPower[™] Ligation PreMix

Linker ligation

QUALITY CONTROL

Cohesive-end Ligation	> 95% of Hind III fragments of Lambda DNA are ligated in 5 minutes when 1 ug of DNA fragments in 20 ul is incubated with AccuPower™ Ligation PreMix at room temperature.
Blunt-end	> 95% of EcoR V fragments of Lambda DNA are ligated in 10 minutes when 1 ug of DNA fragments in 20 ul is incu-
Ligation	bated with AccuPower™ Ligation PreMix at room temperature.

APPLICATIONS

Cloning into vectors

Library construction

• TA cloning

Recirculization of linear DNA

STORAGE CONDITIONS

Storage temperature : -20°C


EXPERIMENTAL DATA



Lane 1:1 ug of λ DNA / Hind III fragment Lane 2, 3: Cohesive-end ligation product with AccuPower™ Ligation PreMix Lane 4, 1 ug of λ DNA / EcoRV fragment


Lane 5, 6: Blunt-end ligation product with AccuPower™ Ligation PreMix

Fig 2. Ligation efficiency test of AccuPower™ Ligation PreMix using various reaction times.



Lane A: 1 ug of λ DNA / Hind III fragment Lane B: 1 ug of λ DNA / EcoRV fragment Lane 1: Ligation incubation for 60 min at RT Lane 2: Ligation incubation for 50 min at RT Lane 3: Ligation incubation for 40 min at RT Lane 4: Ligation incubation for 30 min at RT Lane 5: Ligation incubation for 20 min at RT Lane 6. Ligation incubation for 10 min at RT Lane 7: Ligation incubation for 5 min at RT

Fig 3. Ligation comparison between AccuPower™ Ligation PreMix and other competitor's products



Lane A: Intact \lambda DNA

Lane A: Intact \DNA Lane B: \DNA / Hind III fragment Lane C: \DNA / EcoRV fragment Lane 1, 2: AccuPower^m Ligation PreMix Lane 3: 14 DNA Ligase from company N Lane 4: Quick Ligation Kirform company N Lane 5: LigaFast Rapid DNA Ligation System from company P Lane 6: Ready-To-Go T4 DNA Ligase from company A

ORDERING INFORMATION

PRODUCTS	SIZE	CAT.NO
AccuRapid [™] DNA Ligation Kit	50 ligations	K-7101
AccuRapid [™] DNA Ligation Kit	150 ligations	K-7102
AccuPower™Ligation PreMix	96 ligations	K-7103

CONTACT



BIONEER, INC.

1000 Atlantic Avenue, Alameda, CA 94501-1147 USA Toll Free : (877) 264 - 4300 Fax: (510) 865 - 0350 Email : order.usa@bioneer.com

BIONEER CORPORATION

49-3 Munpyeong - dong Daedeok - gu, Daejeon 306 - 220 South Korea Domestic Phone: 82-42-930-8521 International Phone: 82-42-930-8523

Fax: 82-42-930-8622 Fax: 82-42-930-8688

