

# AccuPower® PCR PreMix

ISO 9001 Certified

AccuPower® PCR PreMix Is covered by Korea Pat. No. 162770 and US Pat. No. 5861251 & 6153412. Headquarters 49-3,Munpyeong-dong, Daedeok-gu,Daejeon306-220, Korea Phone:+82-42-936-8500

Fax:+82-42-930-8600

Seoul Office 1 floor, Donggwan Building 1566-2, Seocho-dong, Seocho-gu, Seoul 137-874 Republic of Korea Phone: +82-2-598-1094 Fax:+82-2-598-1096

Bioneer, Inc.
1000 Atlantic Avenue
Alameda, CA 94501 USA
Toll free: 1-877-264-4300
Fax: 1-510-865-0350
E-mail:order.usa@bioneer.com

Order
Domestic only: 1588-9788
E-mail: order@bioneer.co.kr
URL: www.bioneer.com

AccuPower® PCR PreMix is a new, powerful, and ready-to-use PCR reagent optimized for more accurate PCR amplifications.

## **Advantages**

**Speed** Substantial reduction in PCR setup time.

Stability

AccuPower® PCR PreMix is the powerful technology for convenient and easy to perform DNA amplification. It contains DNA polymerase, dNTPs, a tracking dye and reaction buffer in a premixed format, freeze-dried into a pellet. The patented chemical stabilizer of this product enables to maintain the activity of pre mixture for over a month even when stored at room temperature

(25°C), over 2 years in the freezer.

Reproducibility and Yield Batch manufacturing under strict ISO 9001 quality control conditions guarantees reproducibility. Minimal handling during

reaction set up provides improved accuracy.

Simplicity Each product line contains an application specific enzyme in an easy to re-suspended, lyophilized premix of dNTPs, reaction

buffer, a tracking dye, and a stabilizer. The premix is ready to use, thus offering virtually no reaction set up time.

## Contents

Reaction size Component	20 μl reaction	50 μl reaction
Top DNA polymerase	1 U	2.5 U
Each: dNTP (dATP, dCTP, dGTP, dTTP)	250 μΜ	250 μΜ
Tris-HCI (pH 9.0)	10 mM	10 mM
KCI	30 mM	30 mM
MgCl <sub>2</sub>	1.5 mM	1.5 mM
Stabilizer and tracking dye		

#### **Experimental Protocol**

Add template DNA and primers to AccuPower<sup>®</sup> PCR tube.

Concentration of template DNA and primer

Concentration of template DNA and primer			
	20 μl reaction	50 μl reaction	
Template DNA	5 – 50 ng	10 - 100 ng	
Primer	5 - 10 pmole	10 - 20 pmole	

- Add distilled water to AccuPower® PCR tubes to a total volume of 20 μl or 50 μl.
- Dissolve the lyophilized blue pellet by vortexing, and briefly spin down.
- Add mineral oil. (This step is unnecessary when using a thermal cycler with top heating.)
- · Perform PCR of samples.
- Load samples on agarose gel without adding a loading-dye mixture, and perform electrophoresis.

# **Notice to Purchaser**

This enzyme is specifically optimized for increasing base incorporation rate by inactivating 5'->3' exonuclease activity. Therefore, this is not recommended to use for Real Time PCR using TaqMan® probe.

## Trademark

AccuPower is a trademark of Bioneer Corporation.

## **Ordering Information**

Tube type	Reaction	Cat.No.	Description	
0.2ml Tube -	20 ul	K-2012	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes	
		K-2016	0.2 ml thin-wall 8-strip tubes with attached cap / 480 tubes	
	50 ul	K-2013	0.2 ml thin-wall 8-strip tubes with attached cap / 96 tubes	
		K-2017	0.2 ml thin-wall 8-strip tubes with attached cap / 480 tubes	
0.5ml Tube	50 ul	K-2011	0.5 ml thin-wall tubes with attached cap / 100 tubes	
96-well	10 ul	K-2260-1	thin-wall 96-well flat plate	
		K-2260-2	thin-wall 96-well full-skirted plate	
		K-2260-3	thin-wall 96-well semi-skirted plate	
	20 ul	K-2260-4	thin-wall 96-well flat plate	
		K-2260-5	thin-wall 96-well full-skirted plate	
		K-2260-6	thin-wall 96-well semi-skirted plate	
384-well	5 ul	K-2080-1	thin-wall 384-well full-skirted plate	
	10 ul	K-2080-2	thin-wall 384-well full-skirted plate	
	20 ul	K-2080-3	thin-wall 384-well full-skirted plate	
Master		K-2018	1 ml of 2 X Master mix solution	
Mix	•	K-2018-1	10 ml of 2 X Master mix solution	

A complete product list appears on our web site at www.bioneer.com